

Compare Tuttnauer EZ11-Plus and Midmark M11



Tuttnauer EZ11-Plus

Midmark M11

Features

- Dynamic air-removal technology and a closed door HEPA-filter air-drying system which maintains sterility
- Optional printer to print detailed-history of each cycle and document time, temperature, pressure, etc.
- Automatic download of cycle history to your PC
- Meets ANSI/AAMI STT55 sterilization standards

- Steam-flush pressure pulse air-removal system
- Pre-programmed sterilization cycles
- Alerts users in real-time of any critical system errors

Overall Dimensions (L x W x H)

24.8" x 20.9" x 17.3"
(63 cm x 53 cm x 44 cm)

18.1" x 17.8" x 17.8"
(50 cm (w/plug) x 45.2 cm x 45.2 cm)

Chamber Dimensions

11" x 19.8"
(28 cm x 50.4 cm)

11" x 18"
(28 cm x 45.7 cm)

Slightly larger chamber

Chamber Volume	28.5 Liters Largest capacity in the tabletop market	24.6 Liters
Included Trays	5 Trays Five large- 16.3" x 6.7" x 0.6" (41.3 cm x 16.9 cm x 1.6 cm) 1 Pouch Rack	4 Trays Two large- 15" x 9" x 1.1" (38 cm x 22.9 cm x 2.9 cm) Two small- 15" x 6.6" x 1.1" (38 cm x 16.8 cm x 2.9 cm) Pouch Rack not included
Included Cassettes	Can be used with included tray rack 4 Full Size Standard Cassettes 4 Half Size Standard Cassettes	Requires accessory rack (not included) 4 Full Size Standard Cassettes 4 Small Cassettes
Factory Warranty	2 year parts & labor warranty 10 year chamber warranty	1 year warranty
Shipping Weight	56 kg	59.4 kg